

Excel

!CLICK-&-GO

Lite



Innovation in wheelchair power support systems

Compact powerpack

Simply be placed under the wheelchair

Easily operable



The Excel Click & Go *Lite* is equipped with two drive wheels and programmable microprocessor-driven electronics. The unit can be fitted to all wheelchairs with seat widths between 40 and 55 cm. Easy operating at the push handles. The Click & Go *Lite* is light enough to take anywhere you go, and when you do not need it, the pack is easily lifted of the floor.



Small enough to take everywhere you go.



Always smart: easy to adjust speed forwards and backwards drive.



Easy to fit

‘Smart mobility solutions’

Maximum speed:	4 mph
Range:	8 km
Motor:	12V, 90W-280W
Battery:	12V 20Ah
Charger:	3A/12V

Fits wheelchairs seat width with:	40 to 55 cm
-----------------------------------	-------------



Product complies with crash test EN12184

We have a continuous product development policy and consequently reserve the right to amend design and specifications without prior notice. All errors and omissions excepted. All rights reserved. Copyright 2007 Van Os Medical B.V.

Overall width:	25 cm
Overall length:	44 cm
Overall height:	29½ cm
Drive Wheels:	200x50 mm

Weight (complete):	14 kg
Weight (without batteries):	7,7 kg
Weight (battery):	6,3 kg

The Click & Go comes standard with:

- AccuSet including carrying case
- Battery chargers
- Manual operation cabinet
- Mount



Van Os Medical UK Ltd
42-44 Ousegate, Selby
YO8 4NH

Tel.: +44 (0) 1757 701177
Fax: +44 (0) 1757 706011

sales@vanosmedical.com

www.vanosmedical.com