

Excel

G3 'paediatric'



Lightweight paediatric wheelchair

Equipped with height adjustable push handles

Easy foldable and very compact

Height adjustable push handles.

Foldable, tension adjustable nylon seat and backrest upholstery.

Depth and angle adjustable aluminium brakes.

Lightweight aluminium wheels, fitted with PU tyres.

Lightweight fitted spokes with aluminium hand rims.

Fitted with durable castors 6" x 1¼" (PU).

Backrest height adjustable, 38, 40 and 42 cm.
Foldable aluminium cross bar.

Durable steel side frames, with lightweight aluminium armrests and footrests.

Easy swing inwards / outwards removable footrests.

Tie down points illustrated for transportation.

Adjustable nylon heel loops on the foot plates.



Very easy to fold, very compact and easily transported.



Easy swing inwards/outwards removable footrests.



Flip backwards and removable armrest. Including skirtguard and optional skirtguard with spray guard.



Height adjustable push handles.

'Smart mobility solutions'

Standard specifications

| | |
|------------------|--------------------|
| Seat width: | 30 and 35 cm |
| Seat depth: | 30 and 35 cm |
| Seat height: | 40½ and 43 cm |
| Backrest height: | 38, 40 and 42 cm |
| Total width: | Seat width + 21 cm |
| Total length: | > 83 cm |
| Total height: | > 80 cm |
| Width folded: | 28 cm |
| Armrest height: | 20 cm |

| | |
|--------------------|-----------------------------------|
| Front castor size: | 6" x 1" (other sizes available) |
| Rear wheel (sp): | 22" x 1¾" (other sizes available) |

| | |
|---------------------|---------|
| Total weight: | 13.9 kg |
| Weight (transport): | 9 kg |
| Max. user weight: | 75 kg |

Frame colour available in:



Product complies with crash test ISO 7176-19:2001 and EN12183



Van Os Medical UK Ltd
Corunna House
42-44 Ousegate
Selby
YO8 4NH

Tel.: +44 (0) 1757 701177
Fax: +44 (0) 1757 706011

sales@vanosmedical.com

www.vanosmedical.com

All rights reserved. Copyright 2008© Van Os Medical B.V.®