

Torsion suspension wheelchair

Torsion suspension and independent wheel positioning
Progressive suspension on frontcastors

Backrest with 8 tension straps
Backrest height 47½ and 50 cm (adjustable)

Armrests of comfortable polyurethane,
Long foldable, detachable armrests

Seat cushion, 5 cm with non-skid
layer

Footrests are easily detachable
and height adjustable

Frame with stabilisation

One piece footplate, detachable
which makes it possible to use as to
independent footplates

Height and angle adjustable push handles

Depth and angle adjustable
aluminum brakes

Torsion suspension and
independent wheel positioning
Adjustable suspension

Height adjustable, detachable anti
tip wheels with steppers

Front castors 7" x 1¼" (PU), suspended
and multi position
Angle adjustment on ballheads



One piece footplate, foldable,
can be detached for use as two
separate footplates.



Independent torsion suspension
and independent wheel
positioning for optimal comfort
for the users back.



Suspended multi position
design front castors which
completes this super
suspended wheelchair.

'Smart mobility solutions'


Standaard specifications

Seat width: 40, 45 and 50 cm
Seat depth: 46 cm
Seat height: 49, 51 and 53 cm
Backrest height: 47½ and 50 cm (adjustable)
Total width: Seat width+ 20 cm
Total length: > 101 cm
Total height: > 90 cm
Width folded: 35 cm
Armrest height: 26 cm

Front castor size: 7" x 1¼" (other sizes available)
Rear wheel (sp): 24" x 1" (other sizes available)

Total weight: 19 kg
Weight (transport): 12½ kg
Max. user weight: 130 kg

Frame colour available
in:



Product complies with EN12183

All rights reserved. Copyright 2009 Van Os Medical B.V.®



Van Os Medical UK Ltd
Corunna House
42-44 Ousegate
Selby
YO8 4NH

Tel.: +44 (0)1757 701177
Fax: +44 (0)1757 706011

sales@vanosmedical.com