

Lightweight modular wheelchair

Foldable modular lightweight frame

Also available in pediatric sizes

Half folding backrest can be fitted as standard for UK customers. (For 55 cm seat width, only fixed backrest available).

Backrest equipped with 8 adjustable straps. Backrest height adjustable 42½, 45 and 47½ cm.

Armrests made of durable polyurethane, contoured for extra comfort. Height adjustable, foldable armrests (optional).

Seat cushion, 5 cm deep with non slip base. Adjustable seat upholstery.

Footrests detachable, swing away and swing in. Height adjustable footplates.

Tie down points illustrated for transportation.

Footplate angle adjustable 90 - 105. Adjustable nylon heel loops.

Front castors 7" x 1 1/4" polyurethane. Multi position front fork, angle adjustable front castor, horizontal adjustable.

Height adjustable push handles, also angle adjustable.

Backrest foldable tension bar, offering extra stability and greater comfort for the user.

Depth and angle adjustable braking system with extra comfy levers.

Multi adjustable rear axle plate.

Quick release detachable rear wheel. S/P and transit. Lightweight, fitted with PU tyre.

Anti-tip system fitted as standard. Height adjustable and detachable. With integrated tipper system.

Seat height is adjustable.



Also available in Transit version, fitted with 12" rear wheels. Seat height standard is 52 cm.



With:
- Seat height adjustable.
- Flat free (PU) front and rear wheels.



Supremely comfortable length compensating elevating legrest option. Includes easy adjustable calf support.



Optional high/low adjustable flip backward and detachable armrests.

'Smart mobility solutions'

Standard specifications

Seat width: 35, 37½, 40, 42½, 45, 50 and 55 cm
 Seat depth: 43 and 46 cm
 Seat height: 43, 45, 47, 50, 52 and 54 cm
 Backrest height: 42½, 45 and 47½ cm (adjustable)
 Total width: Seat width + 17 cm
 Total length: > 103 cm
 Total height: > 95 cm
 Width folded: 31 cm
 Armrest height: 26 or 19-35 cm adjustable (optional)

Front castor size: 7" x 1 1/4" (other sizes available)
 Rear wheel (sp): 24" x 1 3/8" (other sizes available)
 Rear wheel (transit): 12½" x 2¼"

Total weight (from): 13.9 kg (depending on the model)
 Weight (transport): 8.4 kg
 Max. user weight: 130 kg (seat width 50 and 55, 150 kg)

Frame colour available in:



Product complies with crash test ISO 7176-19:2001 and EN12183

All rights reserved. Copyright 2008® Van Os Medical B.V.®



Van Os Medical UK Ltd
 Corunna House
 42-44 Ousegate
 Selby
 YO8 4NH

Tel.: +44 (0) 1757 701177
 Fax: +44 (0) 1757 706011

sales@vanosmedical.com