

Lightweight paediatric wheelchair

Aluminium lightweight frame

Seat width and seat depth adjustable 25 to 35 cm



Height and angle adjustable push handles.

Height adjustable and detachable armrests.

Depth and angle adjustable aluminum brakes with extended brake handles.

Wheel position adjustable, with three different positions.

Detachable rear wheels, With easy to use quick-release button aluminum wheels with lightweight aluminum hand rims.

Anti-tip system height adjustable using a quick-release system, also turnable.



Backrest equipped with 5 tension straps.
Adjustable backrest angle (one time adjustment) 90 to 105 degrees.
Seat width and depth adjustable from 25 to 35 cm.
Backrest height adjustable (one time adjustment).
Adjustable aluminum cross brace.

Seat cushion as standard, 5 cm thickness, non-slip base. Tension adjustable seat canvas.

Swing-away footrests, footrest angle 80 degrees.

Tie down points illustrated for transportation.

Adjustable nylon heel loops on the foot plates.

Front castor wheels 6" x 1 1/4" (PU). Standard equipped with stylish front forks.

Growing, changing and learning... The choice of a child's wheelchair is always determined by the size of the child at that time. This will always be the case when you are looking for a suitable short term solution.

But what about the child's needs in the future?

With these thoughts in mind Excel has developed the Excel G5 Paediatric wheelchair that grows with the child.

This allows the child to use the same wheelchair for many years without the need for the costly investment in many different wheelchairs. The G5 Paediatric wheelchair is seat width, seat depth and seat height adjustable, meaning one wheelchair will last the child many years.



With:
- Adjustable nylon heel loops on the foot plates
- Depth and angle adjustable brakes



With:
- Heavy-duty aluminum cross brace
- Anti-tip system
- Height adjustable push handles

Standard specifications

Seat width: 25, 30 and 35 cm
Seat depth: 25, 30 and 35 cm
Seat height: 42½ and 45 cm
Backrest height: 25-35 cm (adjustable)
Total width: Seat width + 22 cm
Total length: > 92 cm
Total height: > 80 cm
Width folded: 34 cm
Armrest height: 16-35 cm adjustable

Front castor size: 6" x 1 1/4" (other sizes available)
Rear wheel (sp): 22" x 1" (other sizes available)

Total weight (complete): 14 kg
Weight (transport): 11 kg
Max. user weight: 75 kg

Frame colour available in:



Product complies with crash test ISO 7176-19:2001



Van Os Medical UK Ltd
Corunna House
42-44 Ousegate
Selby
YO8 4NH

Tel.: +44 (0) 1757 701177
Fax: +44 (0) 1757 706011

sales@vanosmedical.com