

Multi-positional tilt wheelchair

Tilt wheelchair with adjustable backrest

14 adjustment options in one frame

Back coating of Dartex material.
Stepless tilt adjustment, easy operable with the control levers.

Armrests height 20 cm, steplessly adjustable till 28 cm, forward/backward and width adjustable. Equipped with padded side cushions.

Abduction bobbin, height and depth adjustable (optional).

Length compensating elevating legrests, including calf plates.

Depth adjustable brakes.

Wheel base adjustable.

Footplates height, depth and angle adjustable, equipped with adjustable heel loops.

Head support, depth, height, angle and side adjustable.
Head support equipped with adjustable side support.

Height and angle adjustable push bar.

Side support, height and depth adjustable. Also angle adjustable and foldable. Equipped with detachable cover.

Backrest, steplessly height adjustable.

Backrest, steplessly angle adjustable, 0-40 degrees.

Seat depth, standard 40 cm, stepless adjustment to 58 cm

Detachable quick release 24" rear wheel. Lightweight spokes, including high pressure performance tyres.

Equipped with anti-tip wheels, height adjustable and detachable. With integrated tipper system.

Front seat height adjustable, standard 50 cm. (3 positions)



With mechanical tilt and back adjustment for a comfortable seat position.



Length compensating elevating legrest, swing away/swing in, detachable. Footplates depth adjustable and angle adjustable.

'Smart mobility solutions'

Standard specifications

Seat width: 40-45, 45-50 and 50-55 cm
Seat depth: 40 adjustable to 50 cm in one frame
Seat height: 47½, 50 and 52½ cm
Backrest height: 48 - 58 cm (adjustable)
Total width: 68 cm
Total length: > 100 cm
Total height: > 95 cm
Width folded: 31 cm
Armrest height: 20 - 28 cm adjustable

Front castor size: 200 x 30 mm
Rear wheel (sp): 24" x 1"

Total weight (from): 29 kg (depending on the model)
Weight (transport): 16½ kg
Max. user weight: 150 kg

Frame colour available in:



Product complies with crash test ISO 7176-19:2001 and EN12183

All rights reserved. Copyright 2008© Van Os Medical B.V.®



www.vanosmedical.com



Van Os Medical UK Ltd
Corunna House
42-44 Ousegate
Selby
YO8 4NH

Tel.: +44 (0) 1757 701177
Fax: +44 (0) 1757 706011

sales@vanosmedical.com