

## Aluminium transit wheelchair

Affordable wheelchair with a range of high quality features.

Foldable back and seat coating.

Comfortable armrests.

Optional seat cushion, 5 cm with non slip base.

Standard with padded seat, easily folded.

Standard equipped with safety belt.

Quick release detachable and swing away footrests.

Height adjustable footplates, equipped with nylon heel loops.

Equipped with brakes and durable push handles.

Half folding backrest.

Equipped with a tipping lever for climbing pavements.

Standard with maintenance free PU front and 12" PU rear wheels.



Transit wheelchair equipped with a half folding backrest.



Quick release removable and swing away footrests.




Equipped with brakes and durable push handles.

'Smart mobility solutions'

### Standard specifications

Seat width:	45cm
Seat depth:	41cm
Seat height:	50cm
Backrest height:	43cm
Total width:	Seat width + 13cm
Total length:	> 95cm
Total height:	> 93cm
Total weight:	8.9 kg (Depending on model)
Max. User Weight:	100 kg

Available in colour:



We have a continuous product improvement policy and consequently reserve the right to amend design and specifications without prior notice. All errors and omissions excepted.  
All rights reserved. Copyright 2009© Van Os Medical



Van Os Medical UK Ltd

42-44 Ousegate, Selby  
YO8 4NH

Tel.: +44 (0)1757 701177  
Fax: +44 (0)1757 706011

sales@vanosmedical.com