

Excel Navigator 4



All-round suspension mobility scooter

Equipped with an anti-tip system at the rear
4 wheel mobility scooter with extra comfort suspension

Equipped with a rear-view mirror

Easy to use and understand dashboard
With large touch operation buttons
Setup so the user can use both their left and right hand
Theft prevention system using ignition key
Speed adjustable in 6 steps

Indicators fitted, left and right flashing
with beep to inform of your intended
movements

Height-adjustable steering column
Adjustable steering angle, stepless adjustment,
also can be folded down for transportation
purposes

Impact-resistant plastic covers

Equipped with independent
suspension on the front and rear
wheels

Comfortable head support for
optimal use and comfort

Comfortable captains seat with:
An adjustable seat depth
Adjustable backrest
Lumbar adjustment

Folding armrest, width
angle and height adjustable

Removable collapsible
shopping basket
Including crutch holder

Detachable seat four adjustable
positions for seat post
360 degrees rotation of seat

Equipped with a 650 watt motor,
equipped with 4 motor brushes
Motor equipped with
electro-magnetic brake
Equipped with programmable
digital electronics, Dynamic Rhino
110 controller, speed adjustable

Easy to adjust rear suspension



Easy to use waterproof dashboard, with
touch operation.



Equipped with independent suspension on
the front and rear wheels. Also, the rear
suspension is adjustable.



Alloy hubs caps with line profile tyres.



Equipped with rub non-slip mat for
safety whilst getting in and out of the
seat

'Smart mobility solutions'

Standard specifications

Overall width:	68 cm	Total weight:	145 kg (depending on model)
Height rear:	120 cm	Max. User Weight:	220 kg
Height front:	114 cm	Available in colours:	
Overall length:	137 cm (including crutch holder)	Product complies with crash test	EN12184
Seat width:	47 cm		
Seat height:	64-71 cm		
Guardrail height:	20-30 cm		
Total seat width:	58 cm		
Wheelbase:	91 cm		
Ground clearance:	16 cm		
Speed:	Adjustable to 15 km/h (10 mph)*		
Range:	Up to 50 km by 73 ah batteries		

* For UK market scooters are only available with maximum speed 8 mph

We have a continuous product improvement policy and consequently reserve the right to amend design and specifications without prior notice all errors and omissions accepted.

All rights reserved. Copyright 2009 Van Os Medical B.V.®



Van Os Medical UK Ltd
Corunna House
42-44 Ousegate
Selby
YO9 4NH

Tel.: +44 (0) 1757 701177
Fax: +44 (0) 1757 706011

sales@vanosmedical.com

www.vanosmedical.com