

Excel X-power 30



Rear wheel drive power wheelchair

Power wheelchair suitable for inside and outside use
Equipped with 2 motors



Adjustable backrest.

Joystick socket.

2 x 36Ah/12v closed gel batteries (optional).

Front castors: 200 x 50 mm (pneumatic).

Swing away and detachable footrests.

Height adjustable push handles.

Arm pads made of polyurethane.

Angle adjustable backrest, one time adjustment, 0 - 30°.

Suspension on the drive wheels.

Equipped with 2 motors.
Max. speed 5 mph.
Range of 35 km.
(depending on weight and road condition)

Rear wheels 12" aluminium.



Optional lighting, see and be seen.



Height adjustable arm pads.
Adjustable to the wishes of the user.



Battery packs are easy to disassemble and assemble.

'Smart mobility solutions'

Standard specifications


Seat width: 46 cm
Seat depth: 38 - 48 cm
Seat height: 50 cm

Total length: > 102 cm
Total height: > 123 cm
Width folded: 63 cm

Motor: 24V / 320 Watt (2x)
Speed: max. 5 mph
Range: Up to 35 km (depending on weight and road condition)

Total weight (from): 48 kg (incl. batteries depending on the model)
Max. user weight: 125 kg

Frame colour available in:



Product complies with crash test EN12184

All rights reserved. Copyright 2008© Van Os Medical B.V.®



www.vanosmedical.com



Van Os Medical UK Ltd
Corunna House
42-44 Ousegate
Selby
YO8 4NH

Tel.: +44 (0) 1757 701177
Fax: +44 (0) 1757 706011

sales@vanosmedical.com