

Excel X-power 5



Foldable power wheelchair

Compact manoeuvrable and foldable power wheelchair with reinforced steel frame.



Height adjustable and flip backwards armrests.

Foldable backrest.

Easy to disassemble 12V / 36Ah (2x) battery pack, which makes transport easier.

Equipped with 2 motors.
Max. speed 4 mph
Maximum range of 25 km
(depending on weight and road condition)

Standard equipped with height adjustable anti-tip wheels.

Adjustable armrests made of polyurethane

Easy operable joystick. Unit weatherproof.

Cushion made of high quality nylon material. Tension adjustable seat and backrest coating. Cushion 5 cm deep.

Reinforced foldable steel frame. Firm crossbrace, which makes the chair extra stable and rigid.

Front wheels 200 x 50 mm (pneumatic).



Foldable backrest which makes the wheelchair very compact and easy to transport.



Almost completely disassemblable which makes the wheelchair very light and easy to transport.



Battery packs are easy to disassemble and assemble.

'Smart mobility solutions'

Standard specifications

Seat width: 40 and 45 cm
Seat depth: 46 cm
Seat height: 51 cm
Backrest height: Seat width + 16 cm
Total length: > 110 cm
Total height: > 101 cm
Width folded: 33 cm
Armrest height: 23 - 33 cm adjustable
Motor: 24V / 320 Watt (2x)
Speed: 4 mph
Range: Up to 25 km (depending on weight and road condition)

Total weight (from): 64 kg (incl. batteries depending on the model)
Max. user weight: 135 kg

Frame colour available in:



Product complies with crash test EN12184



All rights reserved. Copyright 2008® Van Os Medical B.V.®

www.vanosmedical.com



Van Os Medical UK Ltd
Corunna House
42-44 Ousegate
Selby
YO8 4NH

Tel.: +44 (0) 1757 701177
Fax: +44 (0) 1757 706011

sales@vanosmedical.com