

Lightweight foldable wheelchair

Active wheelchair for permanent use, with foldable frame and backrest.

Equipped with horizontal padded striped back.

Foldable, tension adjustable back and seat upholstery.

Angle adjustable backrest (one time adjustment), 90 to 105 degrees.

Height adjustable and removable armrests.

Seat cushion 5 cm.
Adjustable seat upholstery.

Angle adjustable castor housing.

Depth adjustable footplates.

Footplates angle adjustable 90 to 105 degrees.

Footrests detachable, swing away and swing in.

Standard equipped with nylon heels loops (not pictured).

Front castors, 5" x 1 1/4" (polyurethane).

Height adjustable push handles, also angle adjustable.

Backrest tension bar, foldable and width adjustable.

Depth and angle adjustable braking system with extra comfy levers.

Detachable quick release 24" rear wheel. Lightweight spokes, including high pressure performance tyre.

Standard equipped with anti-tip wheels (not pictured).

Foldable crossbrace.



A unique concept with an easy folding system. This wheelchair can be disassembled to become extremely compact for transportation or storage.



Three seat heights and depths possible in one frame! With a maximum user weight of 150 kg on all seat widths.

'Smart mobility solutions'

Standard specifications

Seat width: 40, 42½, 45 and 50 cm
 Seat depth: 40, 43 and 46 cm
 Seat height: 48, 50 and 52 cm
 Backrest height: 42½, 45 and 47½ cm (adjustable)
 Total width: Seat width + 17 cm
 Total length: > 100 cm
 Total height: > 95 cm
 Width folded: 31 cm
 Armrest height: 19-35 cm adjustable (optional)

Front castor size: 5" x 1 1/4" (other sizes available)
 Rear wheel (sp): 24" x 1" (other sizes available)

Total weight (from): 13.9 kg (depending on the model)
 Weight (transport): 8.2 kg
 Max. user weight: 150 kg

Frame colour available in:



Product complies with crash test EN12183



All rights reserved. Copyright 2008 Van Os Medical B.V.®

www.vanosmedical.com



Van Os Medical UK Ltd
 Corunna House
 42-44 Ousegate
 Selby
 YO8 4NH

Tel.: +44 (0) 1757 701177
 Fax: +44 (0) 1757 706011

sales@vanosmedical.com