

Biodynamic positioning wheelchair

5 seat widths in one frame

With dual mechanism posture adjustment



Wheelchair shown with comfort armrests, adjustable headrests and length compensating elevating legrests

Backrest upholstery - pressure compensating breathable fabric.

Comfortable armpads. Height adjustable and removable.

Positioning adjustable lever for backrest and seat. Can be used by the occupant (optional).

Non slip removable cushion, pressure compensating breathable fabric.

Depth adjustable footplates, equipped with nylon heel loops.

Shock absorber in front fork.

Height adjustable push handles, including contoured grips.

Comfortable armrests, depth adjustable. Optional tray holder available.

Detachable side supports, height adjustable. Removable padded support (optional).

Detachable quick release 24" rear wheels, easy use button. Fitted with super lightweight spokes and high pressure tyres.

Anti-tip system, height adjustable. Fitted with integrated tippers.



Wheelchair shown in comfortable position.



Standard fitted with adjustable tension straps for positioning.



5 sizes in one chair! Width adjustable cross brace 40 cm to 50 cm in just 8 minutes.



Length compensating elevating legrest, with foldaway calf support (optional).

'Smart mobility solutions'

Standard specifications

Seat width: 40, 42½, 45, 47½ and 50 cm
 Seat depth: 43 and 46 cm
 Seat height: 45, 48½, 50, 51½ and 54 cm
 Backrest height: 46½ cm
 Total width: Seat width + 17 cm
 Total length: > 105 cm
 Total height: > 95 cm
 Width folded: 33 cm
 Armrest height: 19-35 cm adjustable

Front castor size: 5" x 1¼" (other sizes available)
 Rear wheel (sp): 24" x 1" (other sizes available)

Total weight (from): 17.9 kg
 Weight (transport): 10.2 kg
 Max. user weight: 130 kg

Frame colour available in:



Product complies with crash test EN12183



Van Os Medical UK Ltd
 Corunna House
 42-44 Ousegate
 Selby
 YO8 4NH

Tel.: +44 (0) 1757 701177
 Fax: +44 (0) 1757 706011

sales@vanosmedical.com